

BioCLIA® 500

Efficient - Precise - Reliable - Compact

System Specifications:

- 72 test results / hour
- 50 parameters CE-certified (see product list)
- Built-in 15 inch PC with touch screen
- 12 on board reagent positions, cooled
- 24 sample positions incl. 1 STAT position
- Barcoded QC and calibrators
- Bi-directional HL7- LIS connection
- Random Access
- Dimension (L×W×H): 500mm × 645mm × 648mm
- Weight: 89kg



Product Specification
BioCLIA® 500
www.bioclia.de