

# BioCLIA® 6500

## Efficient - Precise - Reliable

### System Specifications:

- 100 - 218 test results / hour
- 50 parameters CE-certified (see product list)
- Built-in 15 inch PC with touch screen
- 24 on board reagent positions, cooled
- 96 sample positions incl. 2 STAT positions
- Barcoded QC and calibrators
- Bi-directional HL7- LIS connection
- Random Access
- Dimension (L×W×H): 1280mm × 720mm × 700mm
- Weight: 150kg



Product Specification  
BioCLIA® 6500  
[www.bioclia.de](http://www.bioclia.de)